

Specification	s Mi151
Range pl	-2.00 to 16.00 pH
m¹	±699.9 mV / ±1999 mV
Tem	-20.0 to 120.0°C / -4.0 to 248.0°F
Resolution pl	0.01 pH
m'	***************************************
Tem	0.1°C (0.1°F)
Accuracy pl	
(@ 20°C) m'	=======================================
Tem	±0.4°C / ±0.8°F
Typical EMC pl	
Deviation m	=======================================
Tem	
pH Automatic Calibration	1 or 2 point-calibration, with 7 memorized buffers
Offset Calibration	±1 pH
Slope Calibration	from 80 to 108%
Temperature Compensation	n automatic, from -20.0 to 120.0°C / -4.0 to 248.0°F or manual, without temperature probe
pH Electrode	MA 917B/1 (included)
Temperature Probe	MA 831R (included)
Environment	0 to 50°C / 32 to 122°F; max RH 95%
Input Impedance	10 <sup>12</sup> Ohm
Power supply	12 VDC power adapter (included)
Dimensions	230 x 160 x 95 mm
Weight	0.9 kg



## **Glass Electrode & Temperature Probe**

Choose from our wide selection of pH and ORP electrodes at pages 6 and 39.

## **Custom dual level LCD**

Large and easyto-read Custom dual level LCD Display with simultaneous readings and with user-friendly icons.



## Accessories

MA9001 pH 1.68 buffer solution, 230 mL bottle pH 4.01 buffer solution, 230 mL bottle pH 6.86 buffer solution, 230 mL bottle MA9004 MA9006 MA9007 pH 7.01 buffer solution, 230 mL bottle MA9009 pH 9.18 buffer solution, 230 mL bottle MA9010 pH 10.01 buffer solution, 230 mL bottle Refilling solution for pH electrode, MA9012 230 mL bottle

MA9015 Electrode storage solution, 230 mL MA9016 Electrode cleaning solution, 230 mL



bottle MA9310 12 VDC Adapter, 220 V MA9311 12 VDC Adapter, 110 V MA9315 Electrode Holder

MA917B/1 Glass body, double junction refillable pH electrode

Double junction, gel filled ORP electrode with platinum sensor MA831R Temperature probe

MA921B/1

## **Ordering Information**

Mi151 is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
  MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution
- M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution
- M10016 Sachet Electrode Cleaning Solution
- Graduate Pipet
- 12 VDC Adapter
- · Instruction manual